



Digital titrator cartridge, sodium thiosulfate, stabilized, 0.0250 N

Kat. číslo produktu: 2409301

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

For dissolved oxygen determination by Winkler titration. Digits x 0.1 = mg/L DO in 25 mL. Stabilized sodium thiosulfate formulation

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technické údaje

Názov metódy: Winkler

Platforma: Roztok