



Reagent set, total, nitrogen, colorimetric TKN, air freight

Kat. číslo produktu:

2608600

Nebezpečné

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.



For determination of Total Kjeldahl Nitrogen by the Nessler method. (Requires Digesdahl Digestion Apparatus) Packaged for air freight transportation. 250 tests.

Technické údaje

Názov metódy:	Nessler
Platforma:	Roztok
Počet testov:	250
Shipping Method:	Air Shipment