



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Kat. číslo produktu: 2890200
EUR Cena bez DPH: Kontaktujte nás
Dostupné

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Technické údaje

Jednotky:	100 tests
Najmenší krok:	Steps: 0, 25, 50, 100, 200, 400, 600
Parameter:	Free Chlorine, High Range
Platforma:	Testovacie pružky
Poznámka:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	0 - 600 mg/L Cl ₂