



Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Kat. číslo produktu:

2966700

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Technické údaje

Materiál:	Cellulose
Veľkosť balenia:	25 /pk
Veľkosť póru:	10 µm