



Filter Holder, Magnetic

Kat. číslo produktu:	1352900
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Clear polysulfone funnel with imprinted graduations. Magnetic coupling of funnel to stem for fast filter changes. Used with 47 mm membrane filter for coliform testing and prefiltration of many fluids.

Fits standard manifold systems with 8 rubber stopper connections. Capacity 300 ml (10 oz). Measures 23 x 7 cm (9 x 2.75"). Autoclavable.

Technické údaje

Veľkosť balenia: 1