



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Kat. číslo produktu:	2096250
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Technické údaje

Dĺžka:	150 mm
Kapacita:	Capacity to Neck 50 mL
Priemer:	25 mm
Rozmery (H x D):	25 mm x 150 mm
Špeciálne funkcie:	Tube, 25 x 150 mm
Veľkosť balenia:	10 /pk