



BAG, WHIRL-PAK* 5-1/2X9 PK/10

Kat. číslo produktu:	1437297
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Technické údaje

Kapacita:	720 mL
Veľkosť balenia:	10