



ARSENIC TEST KIT

Kat. číslo produktu: 2822800

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Hach's most economical test kit for determination of arsenic levels.

Simple, economical kit for measuring arsenic levels in the field

Provides fast, reliable results in 20 minutes

Technické údaje

Hmotnosť:	1,13 kg
Model testovacej súpravy:	Test Strip (EZ)‡
Najmenší krok:	Steps: 0, 35, 75, 175, 1,500, 4,000 ppb
Parameter:	Arsenic
Platforma:	Testovacie prúžky
Počet testov:	100
Poznámka:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	0 - 500 ppb As
Rozsah merania (2):	0 - 4000 ppb As
Smallest increment:	Steps: 0, 10, 25, 50, 250, 500 ppb
Štýl kufríka:	D