



## Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Kat. číslo produktu:	6788302
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### Technické údaje

Obsah balenia:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions
Product Kit:	5500 sc Silica Analyser, 2/4 Channels

---

### Obsah balenia

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions