



db CDC304 VODIVOSTNÁ ELEKTRÓDA

Kat. číslo produktu: 945-002

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Vodné a nevodné médiá.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

Technické údaje

Obsah balenia: Includes: Radiometer Analytical CDC304 conductivity cell with Certificate of Conformity.

Obsah balenia

Includes: Radiometer Analytical CDC304 conductivity cell with Certificate of Conformity.