



pHC4006-9 kombinovaná pH elektróda,

Kat. číslo produktu:

E16M333

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Dlhá elektróda s malým priemerom na použitie v sklenených skúmavkách.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

Glass sensors made for long-life and durable use.

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.