



Cell adapter, 12 and 13 mm, for 2100N and 2100AN turbidimeters

Kat. číslo produktu: 3033400

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Allows use of 12 or 13 mm vials with 2100N / AN laboratory turbidimeters, incl. adapter and vial cover

Used with Hach spectrophotometers and portable laboratories. Makes it possible to use other cells in instruments which normally take a 25-mL cell. For use with Models 2100N and AN.

Technické údaje

Hmotnosť:	0,05 kg
Kompatibilita kyviek:	12 mm and 13 mm vials
Obsahuje:	Vial adapter, vial cover, instructions