



Flow-through cell kit, 230 V

Kat. číslo produktu: 4745002

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Low-pressure flow-cell kit allows for a constant optical path. Sample flow controlled by a flow valve module.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells. For rapid, convenient processing of large numbers of samples.

Technické údaje

Certifikáty:	TUV/CE
Hmotnosť:	3,2 kg
Objem vzorky:	120 mL purge volume
Obsahuje:	Flow-cell stand assembly, flow-cell inlet reservoir, reservoir cover, collection drain assembly, flow-cell assembly, interconnecting tubing, flow-cell light shield
Vstupné napätie:	230 V AC
Výstupy:	12V DC at 1 A