



Test kit, phosphate, model PO-23

Kat. číslo produktu:

224902

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Číslo metódy:	Color disc/Ascorbic acid
	for clear samples
Model:	PO-23
Model testovacej súpravy:	PO-23
Parameter:	Phosphate
Platforma:	Farebný disk
Počet testov:	100
Prepravná hmotnosť:	4
Rozsah merania:	0,1-4 mg/L PO ₄
Rozsah merania (2):	1-40 mg/L PO ₄
Smallest increment:	0.1
	1
Štandardný rozsah merania:	0.1 - 4 mg/L PO ₄ ; 1 - 40 mg/L PO ₄
Štýl kufríka:	F