



Millipore SDI Testing Kit

Kat. číslo produktu:

2195900

Nebezpečné

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.



Technické údaje

Obsah balenia:

Kit includes a 47mm diameter filter disc holder, pressure regulator, pressure gauge, PTFE tape, tweezers, membrane filter discs, required valves and fittings, and instruction sheet .

Obsah balenia

Kit includes a 47mm diameter filter disc holder, pressure regulator, pressure gauge, PTFE tape, tweezers, membrane filter discs, required valves and fittings, and instruction sheet .