



## Iron color disc test kit, model IR-18B

**Kat. číslo produktu:**

**146401**

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dostupné



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Technické údaje

|                            |                                    |
|----------------------------|------------------------------------|
| Model:                     | IR-18B                             |
| Model testovacej súpravy:  | IR-18B                             |
| Názov metódy:              | Color disc/<br>1,10 Phenanthroline |
| Parameter:                 | Iron                               |
| Platforma:                 | Farebný disk                       |
| Počet testov:              | 100                                |
| Rozsah merania:            | High Single Range                  |
| Smallest increment:        | 0.25                               |
| Štandardný rozsah merania: | 0.25 - 7 mg/L Fe                   |
| Štýl kufríka:              | B                                  |
| Typ:                       | S jednou stupnicou                 |