



Iron color disc test kit, model IR-18A

Kat. číslo produktu:

146500

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model:	IR-18A
Model testovacej súpravy:	IR-18A
Názov metódy:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	Low Single Range
Smallest increment:	0.02
Štandardný rozsah merania:	0.02 - 1.00 mg/L Fe
Štýl kufríka:	B
Typ:	S jednou stupnicou