



Test kit, molybdate, model MO-2

Kat. číslo produktu:

1419301

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Technické údaje

Model:	MO-2
Model testovacej súpravy:	MO-2
Názov metódy:	Color disc/ Mercaptoacetic acid
Parameter:	Molybdate - As Molybdenum
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	Dual Range
Rozsah merania (2):	0 - 40 mg/L Mo
Smallest increment:	0.2
	1
Štýl kufrika:	F