



Test kit cyanide, model CYN-3

Kat. číslo produktu:

201002

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	CYN-3
Model testovacej súpravy:	CYN-3
Názov metódy:	Color disc/ Pyridine-Pyrazolone
Parameter:	Cyanide
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0 - 0,3 mg/L CN ⁻
Smallest increment:	0.02
Štandardný rozsah merania:	0 - 0.3 mg/L CN ⁻
Štýl kufríka:	R