



Test kit, free chlorine, model CN-70F

Kat. číslo produktu:

1454201

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Quickly and accurately measure free chlorine levels in the field or in the lab

Technické údaje

Model testovacej súpravy:	CN-70F ²
Názov metódy:	Color disc/DPD
Parameter:	Chlorine Free
Platforma:	Farebný disk
Počet testov:	200
Poznámka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	Dual Range
Rozsah merania (2):	0 - 3,4 mg/L Cl ₂
Smallest increment:	0.02
	0.1
Štýl kufríka:	D