



Test Kit, chlorine (free and total), model CN-80

Kat. číslo produktu:
EUR Cena bez DPH:
Dodávka do 2 týždňov

2129000
Kontaktujte nás

Nebezpečné



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	CN-80
Model testovacej súpravy:	CN-80 ²
Názov metódy:	Color Disc/DPD
Parameter:	Chlorine, low range - As free & total Cl ₂
	total Cl ₂
	(0-10 mg/L: T only)
Platforma:	Farebný disk
Počet testov:	100
Prepravná hmotnosť:	2,5
Rozsah merania:	0-0,68 mg/L Cl ₂
Rozsah merania (2):	0 - 3,4 mg/L Cl ₂
Smallest increment:	0.02
	0.1
	0.5
Štandardný rozsah merania:	0 - 0.68 mg/L Cl ₂ ; 0 - 3.4 mg/L Cl ₂ ; 0 - 10 mg/L Cl ₂
Štýl kufríka:	D