



Dissolved oxygen test kit, model OX-2P

Kat. číslo produktu:

146900

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Technické údaje

Model testovacej súpravy:	OX-2P
Názov metódy:	Drop count titration/ Modified Winkler
Parameter:	Dissolved Oxygen - As DO
Platforma:	Drop Count Titration
Počet testov:	100
Rozsah merania:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Smallest increment:	0.2
Štýl kufrika:	D