



Sulfite Test Kit, Model SU-5

Kat. číslo produktu:

148002

EUR Cena bez DPH:

Kontaktujte nás

Dostupné

Nebezpečné



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	SU-5
Názov metódy:	Drop count, titration Iodometric
Parameter:	Sulfite
Platforma:	Drop Count Titration
Počet testov:	100
Rozsah merania:	1.0 - 20 mg/L & 10 - 200 mg/L SO_3^{2-}
Smallest increment:	1
	10
Štýl kufríka:	D