



Test kit, sulfite, model SU-DT

Kat. číslo produktu:

2063300

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	SU-DT
Najmenší krok:	0,4-4,0**
Názov metódy:	Digital Titrator/Iodometric
Parameter:	Sulfite - As SO_3^{2-}
Platforma:	Digitálny titrátor
Počet testov:	100
Rozsah merania:	4 - 800 mg/L SO_3^{2-}
Štýl kufríka:	F