



Test kit, alkalinity, model AL-DT

Kat. číslo produktu:

2063700

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	AL-DT
Najmenší krok:	0.1-10
Názov metódy:	Digital Titrator/Sulfuric acid
Parameter:	Alkalinity - (P) & Total (MO) as CaCO ₃
Platforma:	Digitálny titrátor
Počet testov:	100
Poznámka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Rozsah merania:	10 - 4,000 mg/L CaCO ₃
Štýl kufríka:	F