



Test kit, ozone, model OZ-2

Kat. číslo produktu:

2064400

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

Nebezpečné



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	OZ-2
Model testovacej súpravy:	OZ-2
Názov metódy:	Color disc/DPD
Parameter:	Ozone
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0 - 2,2 mg/L O ₃
Smallest increment:	0.05
Štandardný rozsah merania:	0 - 2.2 mg/L O ₃
Štýl kufríka:	B