



Cylinder, graduated, mixing, glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Kat. číslo produktu: 189641
EUR Cena bez DPH: Kontaktujte nás
Dostupné

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Technické údaje

Certifikácia:	Class B
Kapacita:	50 mL
Materiál skrine:	Glass
Metóda merania:	To Contain
Typ:	S jednou stupnicou, zmiešavací
Veľkosť balenia:	1 each