



## Bottle, BOD, 300 mL, 24/pk

Kat. číslo produktu:

62124

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

---

### Technické údaje

|                   |                    |
|-------------------|--------------------|
| Body Material:    | Glass              |
| Kapacita:         | 300 mL             |
| Materiál:         | Borosilicate Glass |
| Package size UOM: | 24pk               |
| Veľkosť balenia:  | 24 /pk             |