



Be Right™



Beaker, polypropylene, low form, 400 mL, 6/pk

Kat. číslo produktu:

108078

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technické údaje

Jednotky:	6 /pk
Kapacita:	400 mL
Materiál:	Polypropylene, Low Form
Materiál skrine:	Polypropylene
Priemer, vonkajší:	99.1 mm with spout