



Be Right™



Bottle, Storage, Wide Mouth, Write-On, 500 mL, 12/pk

Kat. číslo produktu:

2930003

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

Large dedicated labeling areas help to standardize the way bottle contents are identified. A separate area is provided for name and a larger area allows for detail of contents or safety information and includes a defined space for date information. Bottles are molded of rugged HDPE in 500 and 1000 mL sizes and are available as narrow or wide mouth with leak proof polypropylene cap.

Technické údaje

Body Material:	High Density Plastic
Kapacita:	500 mL
Materiál:	HDPE
Package size UOM:	12pk
Veľkosť balenia:	12 /bag