



pH electrode 5303, gel electrolyte, 3x diaphragms, 0-80 °C, 6 bar, S8 screw head

Kat. číslo produktu:	LZU5303.97.0002
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

pH electrode with gel electrolyte and S8 screw head. 0 - 14 pH, 0 - 80 °C, up to 6 bar. 3x ceramic diaphragms.

Applications: swimming pools, drinking water, gas

scrubbing towers, boiler water, detoxification,

wwtp outlet water, etc.

Technické údaje

Dĺžka:	120 mm
Elektrolyt:	Gel
Káblové spojenie:	S8
Materiál:	Sensor: glass
Membrána:	3 ceramic
Parameter:	pH
Priemer:	12 mm
Reference system:	Encapsulated Ag/AgCl crystals
Rozsah měření:	0 - 14 pH
Rozsah tlaku:	až do 6 bar
Teplotný rozsah:	0 - 80 °C
Záruka:	6 mesiacov