



## **pHC2401-7 Combined pH electrode, Red Rod, ring diaphragm, Type 7**

**Kat. číslo produktu:**

**E16M336**

### **ZASTARANÁ POLOŽKA**

**Táto položka už nie je dostupná.**

Navrhované náhrady:

- E16M400
- 945-252

### **Long-term accuracy and reproducibility with Red Rod technology**

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: annular ring. With plug Type 7.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

### **Reproducible measuring conditions ensured by a special bayonet head**

### **Wide temperature range from -10 to 100°C**

### **Fast and accurate response**

### **Stable signal**

### **Trouble-free liquid junction**