



MultiTest kit case (180 x 356 x 188 mm), blue polypropylene

Kat. číslo produktu:	4660500
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Technické údaje

Rozmery (V x Š x H):	180 mm x 356 mm x 188 mm
Štýl kufrika:	G
Výška:	7 in