



Titration Application Pack Generic Chemistry: ORP I>0

Kat. číslo produktu: AP0020.AT1112

EUR Cena bez DPH: Kontaktujte nás

Dodávka do 2 týždňov

Technické údaje

Obsah balenia:	Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle
Parameter:	ORP I>0
Záruka:	24 mesiacov

Obsah balenia

Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle