



Ammonium/Potassium Combination Standard Solution: 5.000 mg/L NH₄-N, 20.000 mg/L K, 1000 mL

Kat. číslo produktu:	LCW892
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Ammonium/Potassium standard for validation purpose of ammonium process probes.

Technické údaje

Koncentrácia:	5000 mg/L NH ₄ -N, 20.000 mg/L K
Parameter:	Ammonium
Prístroj:	Amtax sc