



Reagent powder pillows, ChromaVer[®] 3, chromium, , 5 or 10 mL, pk/1000

Kat. číslo produktu:

1271028

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov



For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Technické údaje

Číslo metódy:	8023
Názov metódy:	1,5-Diphenylcarbohydrazide
Parameter:	Chromium
Platforma:	Práškové reagentie
Počet testov:	1000
Prístroj:	All except DR/3000, DR/2000, DR 900
Rozsah merania:	0,010 až 0,700 mg/L Cr ⁶⁺
Soulad s EPA:	Áno