



Flask, volumetric, polypropylene, 250 mL

Kat. číslo produktu: 1406046

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Technické údaje

Certifikácia:	Not Applicable
Jednotky:	1 each
Kapacita:	250 mL
Materiál skrine:	Polypropylene