



## Flask, volumetric class A, glass, 50 mL

|                           |                 |
|---------------------------|-----------------|
| Kat. číslo produktu:      | 1457441         |
| EUR Cena bez DPH:         | Kontaktujte nás |
| Dátum odoslania neuvedené |                 |

Calibrated “To Contain” and manufactured to ASTM E438, “Glass in Laboratory Apparatus” requirements. Stopper not included—see page 375.

---

### Technické údaje

|                  |               |
|------------------|---------------|
| Certifikácia:    | No            |
| Jednotky:        | 1 each        |
| Kapacita:        | 50 mL         |
| Materiál skrine: | Class A Glass |