



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Kat. číslo produktu: 1472459
EUR Cena bez DPH: Kontaktujte nás
Dodávka do 2 týždňov

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Technické údaje

Body Material:	Plastic
Kapacita:	500 mL
Materiál:	PP
Package size UOM:	12pk
Veľkosť balenia:	12 /pk