



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 1 L, 6/pk

Kat. číslo produktu:	1472463
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Technické údaje

Body Material:	Plastic
Kapacita:	1000 mL
Materiál:	PP
Package size UOM:	6pk
Veľkosť balenia:	6 /pk