



Hydrazine color disc test kit, model HY-2

Kat. číslo produktu:

184900

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model:	HY-2
Model testovacej súpravy:	HY-2
Názov metódy:	Color disc/ p-Dimethylamino- benzaldehyde
Parameter:	Detergents
Platforma:	Farebný disk
Počet testov:	300
Rozsah merania:	0 - 1 mg/L N_2H_4
Smallest increment:	0
Štandardný rozsah merania:	0 - 1 mg/L N_2H_4
Štýl kufríka:	D