



Cyanuric acid test kit, model CY-3

Kat. číslo produktu: 185102
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Aplikácie:	Pool and Spa
Model:	CY-3
Model testovacej súpravy:	CY-3
Názov metódy:	Turbidimetry
Parameter:	Cyanuric Acid
Počet testov:	100
Rozsah merania:	20 - 100 mg/L
Smallest increment:	20
Štýl kufríka:	D