



Cylinder, graduated, mixing, glass, 25 mL ± 0.3 mL, 0.5 mL divisions, glass stopper #13

Kat. číslo produktu: 189640
EUR Cena bez DPH: Kontaktujte nás
Dodávka do 2 týždňov

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Technické údaje

Certifikácia:	Class B
Kapacita:	25 mL
Materiál skrine:	Glass
Metóda merania:	To Contain
Typ:	S jednou stupnicou, zmiešavací
Veľkosť balenia:	1 each