



Bottle, BOD, black, 300 mL

Kat. číslo produktu:	2086047
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Technické údaje

Body Material:	Glass
Kapacita:	300 mL
Materiál:	Borosilicate Glass
Package size UOM:	1pk
Veľkosť balenia:	1 each