



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 500 mL

Kat. číslo produktu: 2089649

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Navrhované náhrady:

- 2089679

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Technické údaje

Jednotky:	1 each
Kapacita:	500 mL
Materiál:	Glass Screw Cap