



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Kat. číslo produktu:

2089676

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Technické údaje

Jednotky:	6/pk
Kapacita:	250 mL
Materiál:	Glass Screw Cap