



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 250 mL, 4/pk

Kat. číslo produktu:	2089876
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Technické údaje

Jednotky:	4/pk
Kapacita:	250 mL
Materiál:	Polymethylpentene
Typ:	Erlenmeyer