



Pipeta, sérologická, sterilná, polystyrén, 5 mL, modrá, 200/bal.

Kat. číslo produktu:	2092667
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Sterilná jednorazová plastová pipeta. Jednotlivo zabalené. Farebné značenie. Kalibrované „na vyliatie“ (TD). Zhoda s normou ASTM E934.

Graduations are calibrated for accurate dispensing to within $\pm 2\%$

Pipets meet BST/TSE requirements; latex-free

DNase/RNase-free, non-pyrogenic, non-cytotoxic, non-haemolytic

Technické údaje

Kapacita:	5 mL
Negative Graduations:	2.0 mL
Stupnice:	0.10 mL
Veľkosť balenia:	200 /bal.