



Cylinder, graduated, polymethylpentene, 25 mL ± 0.3 mL, 0.5 mL divisions

Kat. číslo produktu: 217240

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Technické údaje

Certifikácia:	Class B
Kapacita:	25 mL
Materiál skrine:	Polymethylpentene
Metóda merania:	To Deliver/To Contain
Typ:	S jednou stupnicou
Veľkosť balenia:	1 each