



Reagent and apparatus set, water in oil

Kat. číslo produktu:

2237300

Nebezpečné

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.



Model WO-1 Test Kit

The Model WO-1 test kit provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance. The kit has proven particularly useful on board ships. Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated.

Measures low levels of water in oil

Provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance

Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated

Technické údaje

Model testovacej súpravy:	WO-1
Názov metódy:	Volumetric displacement/ Calcium hydride
Parameter:	Water-in-oil - As %
Počet testov:	25
Prepravná hmotnosť:	3,5
Rozsah merania:	0,05-1 %
Rozsah merania (2):	0-10 %
Smallest increment:	0.05%
	0.20%
Štýl kufríka:	no case